




PHC1 Antibody

Product Code	CSB-PA017891HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78364
Immunogen	Recombinant Human Polyhomeotic-like protein 1 protein (142-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyhomeotic-like protein 1 (hPH1) (Early development regulatory protein 1), PHC1, EDR1 PH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PHC1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA017891HA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.