



PHC1 Antibody, Biotin conjugated

Product Code	CSB-PA017891HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78364
Immunogen	Recombinant Human Polyhomeotic-like protein 1 protein (142-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyhomeotic-like protein 1 (hPH1) (Early development regulatory protein 1), PHC1, EDR1 PH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PHC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.