





PMAIP1 Antibody

Product Code	CSB-PA018229LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13794
Immunogen	Recombinant Human Phorbol-12-myristate-13-acetate-induced protein 1 protein (1-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phorbol-12-myristate-13-acetate-induced protein 1 (PMA-induced protein 1) (Immediate-early-response protein APR) (Protein Noxa), PMAIP1, NOXA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PMAIP1


Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA018229LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA018229LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA018229LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.