



# FAU Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA01825C0Rb
<b>Abbreviation</b>	40S ribosomal protein S30
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P62861
<b>Immunogen</b>	Recombinant Human 40S ribosomal protein S30 protein (1-59AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	40S ribosomal protein S30, Small ribosomal subunit protein eS30, FAU
<b>Species</b>	Human
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	FAU