



FAU Antibody, Biotin conjugated

Product Code	CSB-PA01825D0Rb
Abbreviation	40S ribosomal protein S30
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62861
Immunogen	Recombinant Human 40S ribosomal protein S30 protein (1-59AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S30, Small ribosomal subunit protein eS30, FAU
Species	Human
Research Area	Epigenetics and Nuclear Signaling
Target Names	FAU