



PRMT3 Antibody, Biotin conjugated

Product Code	CSB-PA018733LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60678
Immunogen	Recombinant Human Protein arginine N-methyltransferase 3 protein (311-529AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein arginine N-methyltransferase 3 (EC 2.1.1.-) (Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 3), PRMT3, HRMT1L3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRMT3