



PRMT3 Antibody, Biotin conjugated

Product Code	CSB-PA018733LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60678
Immunogen	Recombinant Human Protein arginine N-methyltransferase 3 protein (311-529AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein arginine N-methyltransferase 3 (EC 2.1.1.-) (Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 3), PRMT3, HRMT1L3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRMT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.