



PROM1 Antibody, HRP conjugated

Product Code	CSB-PA018751LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43490
Immunogen	Recombinant Human Prominin-1 protein (508-792AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prominin-1 (Antigen AC133) (Prominin-like protein 1) (CD antigen CD133), PROM1, PROML1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PROM1