



PROM1 Antibody, HRP conjugated

Product Code	CSB-PA018751LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43490
Immunogen	Recombinant Human Prominin-1 protein (508-792AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prominin-1 (Antigen AC133) (Prominin-like protein 1) (CD antigen CD133), PROM1, PROML1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PROM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.