



PSEN1 Antibody, FITC conjugated

Product Code	CSB-PA018846EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49768
Immunogen	Recombinant Human Presenilin-1 protein (4-218AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Presenilin-1 (PS-1) (EC 3.4.23.-) (Protein S182) [Cleaved into: Presenilin-1 NTF subunit; Presenilin-1 CTF subunit; Presenilin-1 CTF12 (PS1-CTF12)], PSEN1, AD3 PS1 PSNL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PSEN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.