



RAD51 Antibody, Biotin conjugated

Product Code	CSB-PA019268LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06609
Immunogen	Recombinant Human DNA repair protein RAD51 homolog 1 protein (1-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA repair protein RAD51 homolog 1 (HsRAD51) (hRAD51) (RAD51 homolog A), RAD51, RAD51A RECA
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RAD51