



RBPJ Antibody, Biotin conjugated

Product Code	CSB-PA019486LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06330
Immunogen	Recombinant Human Recombining binding protein suppressor of hairless protein (243-485AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Recombining binding protein suppressor of hairless (CBF-1) (J kappa-recombination signal-binding protein) (RBP-J kappa) (RBP-J) (RBP-JK) (Renal carcinoma antigen NY-REN-30), RBPJ, IGKJRB IGKJRB1 RBPJK RBPSUH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RBPJ