



RHD Antibody, Biotin conjugated

Product Code	CSB-PA019677LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02161
Immunogen	Recombinant Human Blood group Rh(D) polypeptide protein (32-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Blood group Rh(D) polypeptide (RHXIII) (Rh polypeptide 2) (RhPII) (Rhesus D antigen) (CD antigen CD240D), RHD
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RHD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.