



RNF4 Antibody, Biotin conjugated

Product Code	CSB-PA019889LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78317
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF4 protein (1-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF4 (EC 2.3.2.27) (RING finger protein 4) (RING-type E3 ubiquitin transferase RNF4) (Small nuclear ring finger protein) (Protein SNURF), RNF4, SNURF
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RNF4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.