



LTBP4 Antibody, Biotin conjugated

Product Code	CSB-PA02049D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2S1
Immunogen	Recombinant Human Latent-transforming growth factor beta-binding protein 4 protein (977-1180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Latent-transforming growth factor beta-binding protein 4 (LTBP-4), LTBP4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LTBP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.