



RTEL1 Antibody, HRP conjugated

Product Code	CSB-PA020563LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ71
Immunogen	Recombinant Human Regulator of telomere elongation helicase 1 protein (975-1219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of telomere elongation helicase 1 (EC 3.6.4.12) (Novel helicase-like), RTEL1, C20orf41 KIAA1088 NHL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RTEL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.