



RTEL1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA020563LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9NZ71 |
| Immunogen | Recombinant Human Regulator of telomere elongation helicase 1 protein (975-1219AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Regulator of telomere elongation helicase 1 (EC 3.6.4.12) (Novel helicase-like), RTEL1, C20orf41 KIAA1088 NHL |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | RTEL1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |