



S100G Antibody, HRP conjugated

Product Code	CSB-PA020644LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29377
Immunogen	Recombinant Human Protein S100-G protein (1-66AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-G (Calbindin-D9k) (S100 calcium-binding protein G) (Vitamin D-dependent calcium-binding protein, intestinal) (CABP), S100G, CABP9K CALB3 S100D
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	S100G