



S100G Antibody, Biotin conjugated

Product Code	CSB-PA020644LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29377
Immunogen	Recombinant Human Protein S100-G protein (1-66AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-G (Calbindin-D9k) (S100 calcium-binding protein G) (Vitamin D-dependent calcium-binding protein, intestinal) (CABP), S100G, CABP9K CALB3 S100D
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	S100G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.