



S100P Antibody, HRP conjugated

Product Code	CSB-PA020645LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25815
Immunogen	Recombinant Human Protein S100-P protein (17-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-P (Migration-inducing gene 9 protein) (MIG9) (Protein S100-E) (S100 calcium-binding protein P), S100P, S100E
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	S100P
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.