



SETD1A Antibody, Biotin conjugated

Product Code	CSB-PA021101HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15047
Immunogen	Recombinant Human Histone-lysine N-methyltransferase SETD1A protein (188-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone-lysine N-methyltransferase SETD1A (EC 2.1.1.43) (Lysine N-methyltransferase 2F) (SET domain-containing protein 1A) (hSET1A) (Set1/Ash2 histone methyltransferase complex subunit SET1), SETD1A, KIAA0339 KMT2F SET1 SET1A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SETD1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.