



SIRPB1 Antibody, Biotin conjugated

Product Code	CSB-PA021335LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00241
Immunogen	Recombinant Human Signal-regulatory protein beta-1 protein (30-371AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal-regulatory protein beta-1 (SIRP-beta-1) (CD172 antigen-like family member B) (CD antigen CD172b), SIRPB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SIRPB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.