



SLC2A3 Antibody, HRP conjugated

Product Code	CSB-PA021553LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11169
Immunogen	Recombinant Human Solute carrier family 2, facilitated glucose transporter member 3 protein (205-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 2, facilitated glucose transporter member 3 (Glucose transporter type 3, brain) (GLUT-3), SLC2A3, GLUT3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC2A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.