



SMARCA4 Antibody, FITC conjugated

Product Code	CSB-PA021801LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51532
Immunogen	Recombinant Human Transcription activator BRG1 protein (808-1092AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription activator BRG1 (EC 3.6.4.-) (ATP-dependent helicase SMARCA4) (BRG1-associated factor 190A) (BAF190A) (Mitotic growth and transcription activator) (Protein BRG-1) (Protein brahma homolog 1) (SNF2-beta) (SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 4), SMARCA4, BAF190A BRG1 SNF2B SNF2L4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SMARCA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.