



SOX14 Antibody, Biotin conjugated

Product Code	CSB-PA022422LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95416
Immunogen	Recombinant Human Transcription factor SOX-14 protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor SOX-14 (Protein SOX-28), SOX14, SOX28
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SOX14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.