



SS18L1 Antibody, FITC conjugated

Product Code	CSB-PA022701LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75177
Immunogen	Recombinant Human Calcium-responsive transactivator protein (47-396AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-responsive transactivator (SS18-like protein 1) (SYT homolog 1), SS18L1, CREST KIAA0693
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SS18L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.