



SS18L1 Antibody, Biotin conjugated

Product Code	CSB-PA022701LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75177
Immunogen	Recombinant Human Calcium-responsive transactivator protein (47-396AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-responsive transactivator (SS18-like protein 1) (SYT homolog 1), SS18L1, CREST KIAA0693
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SS18L1