



SSTR2 Antibody, FITC conjugated

Product Code	CSB-PA022725LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30874
Immunogen	Recombinant Human Somatostatin receptor type 2 protein (304-369AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Somatostatin receptor type 2 (SS-2-R) (SS2-R) (SS2R) (SRIF-1), SSTR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SSTR2