




TAF6 Antibody

Product Code	CSB-PA023095LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49848
Immunogen	Recombinant Human Transcription initiation factor TFIID subunit 6 protein (627-677AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription initiation factor TFIID subunit 6 (RNA polymerase II TBP-associated factor subunit E) (Transcription initiation factor TFIID 70 kDa subunit) (TAF(II)70) (TAFII-70) (TAFII70) (Transcription initiation factor TFIID 80 kDa subunit) (TAF(II)80) (TAFII-80) (TAFII80), TAF6, TAF2E TAFII70
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TAF6

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA023095LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.