





TAS2R38 Antibody

Product Code	CSB-PA023155LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59533
Immunogen	Recombinant Human Taste receptor type 2 member 38 protein (212-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 2 member 38 (T2R38) (PTC bitter taste receptor) (Taste receptor type 2 member 61) (T2R61), TAS2R38, PTC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS2R38

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA023155LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA023155LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.