



# TCL1A Antibody, HRP conjugated

|                     |                                                                                                         |
|---------------------|---------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA023313LB01HU                                                                                      |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                           |
| Uniprot No.         | P56279                                                                                                  |
| Immunogen           | Recombinant Human T-cell leukemia/lymphoma protein 1A protein (1-90AA)                                  |
| Raised In           | Rabbit                                                                                                  |
| Species Reactivity  | Human                                                                                                   |
| Tested Applications | ELISA                                                                                                   |
| Form                | Liquid                                                                                                  |
| Conjugate           | HRP                                                                                                     |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                        |
| Purification Method | >95%, Protein G purified                                                                                |
| Isotype             | IgG                                                                                                     |
| Clonality           | Polyclonal                                                                                              |
| Alias               | T-cell leukemia/lymphoma protein 1A (Oncogene TCL-1) (Oncogene TCL1)<br>(Protein p14 TCL1), TCL1A, TCL1 |
| Immunogen Species   | Homo sapiens (Human)                                                                                    |
| Research Area       | Epigenetics and Nuclear Signaling                                                                       |
| Target Names        | TCL1A                                                                                                   |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                             |