



TLR1 Antibody, FITC conjugated

Product Code	CSB-PA023599LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15399
Immunogen	Recombinant Human Toll-like receptor 1 protein (25-580AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 1 (Toll/interleukin-1 receptor-like protein) (TIL) (CD antigen CD281), TLR1, KIAA0012
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TLR1