



TLR1 Antibody, Biotin conjugated

Product Code	CSB-PA023599LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15399
Immunogen	Recombinant Human Toll-like receptor 1 protein (25-580AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 1 (Toll/interleukin-1 receptor-like protein) (TIL) (CD antigen CD281), TLR1, KIAA0012
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TLR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.