




Tph1 Antibody

Product Code	CSB-PA024100HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17532
Immunogen	Recombinant Mouse Tryptophan 5-hydroxylase 1 protein (1-447AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptophan 5-hydroxylase 1 (EC 1.14.16.4) (Tryptophan 5-monoxygenase 1), Tph1, Tph
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Tph1

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA024100HA01MO at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.