



Tph1 Antibody, Biotin conjugated

Product Code	CSB-PA024100HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17532
Immunogen	Recombinant Mouse Tryptophan 5-hydroxylase 1 protein (1-447AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptophan 5-hydroxylase 1 (EC 1.14.16.4) (Tryptophan 5-monoxygenase 1), Tph1, Tph
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Tph1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.