



TRAF2 Antibody, Biotin conjugated

Product Code	CSB-PA024147LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12933
Immunogen	Recombinant Human TNF receptor-associated factor 2 protein (147-281AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TNF receptor-associated factor 2 (EC 2.3.2.27) (E3 ubiquitin-protein ligase TRAF2) (RING-type E3 ubiquitin transferase TRAF2) (Tumor necrosis factor type 2 receptor-associated protein 3), TRAF2, TRAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRAF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.