



TUBA1C Antibody, HRP conjugated

Product Code	CSB-PA025308HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQE3
Immunogen	Recombinant Human Tubulin alpha-1C chain protein (29-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin alpha-1C chain (Alpha-tubulin 6) (Tubulin alpha-6 chain) [Cleaved into: Detyrosinated tubulin alpha-1C chain], TUBA1C, TUBA6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TUBA1C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.