



TUBB4A Antibody, Biotin conjugated

Product Code	CSB-PA025324LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04350
Immunogen	Recombinant Human Tubulin beta-4A chain protein (1-444AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin beta-4A chain (Tubulin 5 beta) (Tubulin beta-4 chain), TUBB4A, TUBB4 TUBB5
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TUBB4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.