



USP2 Antibody, Biotin conjugated

Product Code	CSB-PA025710LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75604
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 2 protein (9-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 2 (EC 3.4.19.12) (41 kDa ubiquitin-specific protease) (Deubiquitinating enzyme 2) (Ubiquitin thioesterase 2) (Ubiquitin-specific-processing protease 2), USP2, UBP41
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.