



VDAC2 Antibody, FITC conjugated

Product Code	CSB-PA025824LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45880
Immunogen	Recombinant Human Voltage-dependent anion-selective channel protein 2 protein (2-36AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent anion-selective channel protein 2 (VDAC-2) (hVDAC2) (Outer mitochondrial membrane protein porin 2), VDAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VDAC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.