




WRB Antibody

Product Code	CSB-PA026148LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00258
Immunogen	Recombinant Human Tail-anchored protein insertion receptor WRB protein (38-99AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:2000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tail-anchored protein insertion receptor WRB (Congenital heart disease 5 protein) (Tryptophan-rich basic protein) (WRB), WRB, CHD5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	WRB

Image



Western Blot

Positive WB detected in: Mouse spleen tissue


All lanes: WRB antibody at 3.4µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 20, 17 kDa


Observed band size: 20 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA026148LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA026148LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA026148LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.