



# WRB Antibody, Biotin conjugated

Product Code	CSB-PA026148LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00258
Immunogen	Recombinant Human Tail-anchored protein insertion receptor WRB protein (38-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tail-anchored protein insertion receptor WRB (Congenital heart disease 5 protein) (Tryptophan-rich basic protein) (WRB), WRB, CHD5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	WRB