



ZNF2 Antibody, FITC conjugated

Product Code	CSB-PA026572LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSG1
Immunogen	Recombinant Human Zinc finger protein 2 protein (134-289AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 2 (Zinc finger protein 2.2) (Zinc finger protein 661), ZNF2, ZNF661
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.