



RAC1 Antibody, Biotin conjugated

Product Code	CSB-PA04329D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63000
Immunogen	Recombinant Human Ras-related C3 botulinum toxin substrate 1 protein (11-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related C3 botulinum toxin substrate 1 (Cell migration-inducing gene 5 protein) (Ras-like protein TC25) (p21-Rac1), RAC1, TC25
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.