




SUMO2 Antibody

Product Code	CSB-PA05009A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61956
Immunogen	Recombinant Human Small ubiquitin-related modifier 2 protein (1-84AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small ubiquitin-related modifier 2 (SUMO-2) (HSMT3) (SMT3 homolog 2) (SUMO-3) (Sentrin-2) (Ubiquitin-like protein SMT3B) (Smt3B), SUMO2, SMT3B SMT3H2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SUMO2

Image



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA05009A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.