



SUMO2 Antibody, Biotin conjugated

Product Code	CSB-PA05009D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61956
Immunogen	Recombinant Human Small ubiquitin-related modifier 2 protein (1-84AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small ubiquitin-related modifier 2 (SUMO-2) (HSMT3) (SMT3 homolog 2) (SUMO-3) (Sentrin-2) (Ubiquitin-like protein SMT3B) (Smt3B), SUMO2, SMT3B SMT3H2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SUMO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.