





AR Antibody

Product Code	CSB-PA05979A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Recombinant Human Androgen receptor protein (660-919AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Androgen receptor (Dihydrotestosterone receptor) (Nuclear receptor subfamily 3 group C member 4), AR, DHTR NR3C4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AR

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA05979A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA05979A0Rb at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.