





AR Antibody

Product Code	CSB-PA05979A0Rb
Abbreviation	Androgen receptor
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Recombinant Human Androgen receptor protein (660-919AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Relevance	Steroid hormone receptors are ligand-activated transcription factors that regulate eukaryotic gene expression and affect cellular proliferation and differentiation in target tissues. Transcription factor activity is modulated by bound coactivator and corepressor proteins. Transcription activation is down-regulated by NR0B2. Activated, but not phosphorylated, by HIPK3 and ZIPK/DAPK3.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Androgen receptor (Dihydrotestosterone receptor) (Nuclear receptor subfamily 3 group C member 4), AR, DHTR NR3C4
Species	Human
Research Area	Signal Transduction
Target Names	AR

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA05979A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA05979A0Rb at dilution of 1:100