



AR Antibody, HRP conjugated

Product Code	CSB-PA05979B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10275
Immunogen	Recombinant Human Androgen receptor protein (660-919AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Androgen receptor (Dihydrotestosterone receptor) (Nuclear receptor subfamily 3 group C member 4), AR, DHTR NR3C4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.