



IL17A Antibody, Biotin conjugated

Product Code	CSB-PA06659D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16552
Immunogen	Recombinant Human Interleukin-17A protein (26-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-17A (IL-17) (IL-17A) (Cytotoxic T-lymphocyte-associated antigen 8) (CTLA-8), IL17A, CTLA8 IL17
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL17A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.